

## Quotation Ref: QUO-10219-R7C8X6

18 January 2021

Page 1

Animal and Plant Health Agency Woodham Lane New Haw Addlestone Surrey KT15 3NB

Dear Dan,

Please find attached our revised quotation for the various IVC equipment requested. I have grouped this into the following sections:

Table 1 – GM500 items

**Table 2** – 1284 Lids

Table 3 – Autoclavable connections & Logistics

Table 4 – DGM racks, connections, AHU and Change Station

I have changed the 2 x complete DGM80 racks, to 2 x complete discounted DGM70 racks. This includes 140 x Polysulfone cage bases with silicone seal, 140 x stainless steel half wire bar lids, 140 x Polysulfone cage tops, 140 x 260ml Polysulfone bottle with silicone seal, 140 x 25mm bottle caps and 140 x label holders. Cost for these racks would have been over £15,000, each, however, I am pleased to be able to offer this to you for £13,000 total. All racking is still in its original packaging and has only been stored in a warehouse and not been to another site.

All other equipment in this table group remains the same from the original quote.

I will contact you later this week for your feedback.

Many thanks

More details can be found at <u>www.tecniplast.it/uk/</u>

c.c Sam Watson

Tecniplast UK Ltd BCM Box 3058 London, WC1N 3XX

Telephone 0345 050 4556 Fax 0345 050 4557 Email info@tecniplastuk.com





Registered Office: 27 Old Gloucester Street, London WC1N 3AX Registered in England No. 3243587



Page | 2

# The Tecniplast GM500 Caging System

# **Key Benefits Include**

- The chosen IVC for breeding, experimental and quarantine areas throughout the world
- Safe for your staff silicone gasket where 2 points meet
- Maximum Flexibility Select either positive or negative pressure regime to suit your health status
- Guaranteed protection for you and your animals H14 graded HEPA filtration
- Easy to maintain through innovative rack plenum design
- **Optimum use of holding room space** with the smallest footprint of IVC racks available in line with the latest EU guidelines
- Zero cross contamination across all cages on each IVC rack

# Caging options available

- In-cage enrichment
  - o GYM500 Mouse wheel
  - Mouse House
  - o Mouse Mezzanine
  - o Mouse Loft
  - Mouse Balcony
- Multi coloured labels and cage components to distinguish between multiple barriers
- GM500 Plus base suitable for 2017 EU Guidelines for Inbred Trio breeding or monogamous pair with additional female
- GM500 in-cage filter to reduce rack maintenance. Zero consumables and ideal if your facility does not have a rack washer

# **Technical Specification**

• (All rack sizes below are also available as 9 high racks but with the same footprint, please ask for more details):

Single Sided		Double Sided			
DGM60	DGM70	DGM80	DGM120	DGM140	DGM160
1314(w) x 500 (d) x 1998 (*2076) (h)	1528(w) x 500 (d) x 1998 (*2076) (h)	1742 (w) x 500 (d) x 1998 (*2076) (h)	1314 (w) x 888 (d) x 1998 (*2081) (h)	1528 (w) x 888 (d) x 1998 (*2081) (h)	1742 (w) x 888 (d) x 1998 (*2081) (h)
6(w) x10(h)	7(w) x10(h)	8(w) x10(h)	6(w) x10(h)	7(w) x10(h)	8(w) x10(h)
H-Temp, (Polysulfone) 3 years Warranty X-Temp, (Polyphenylsulfone) 8 years Warranty					
Positive and Negative					
	1314(w) x 500 (d) x 1998 (*2076) (h)	DGM60 DGM70   1314(w) x 500 (d) x 1998 (*2076) (h) 1528(w) x 500 (d) x 1998 (*2076) (h)   6(w) x10(h) 7(w) x10(h)   6(w) x10(h) 7(w) x10(h)	DGM60 DGM70 DGM80   1314(w) x 500 (d) x 1998 (*2076) (h) 1528(w) x 500 (d) x 1998 (*2076) (h) 1742 (w) x 500 (d) x 1998 (*2076) (h)   6(w) x10(h) 7(w) x10(h) 8(w) x10(h)   6(w) x10(h) 7(w) x10(h) 8(w) x10(h)   H-Temp, (Polysulfone X-Temp, (Polyphenylsulf 1000	DGM60 DGM70 DGM80 DGM120   1314(w) x 500 (d) x 1998 (*2076) (h) 1528(w) x 500 (d) x 1998 (*2076) (h) 1742 (w) x 500 (d) x 1998 (*2076) (h) 1314 (w) x 888 (d) x 1998 (*2076) (h)   6(w) x10(h) 7(w) x10(h) 8(w) x10(h) 6(w) x10(h)   H-Temp, (Polysulfone) 3 years Warrow X-Temp, (Polyphenylsulfone) 8 years Warrow X-Temp, (Polyphenylsulfone) 8 years Warrow	DGM60 DGM70 DGM80 DGM120 DGM140   1314(w) x 500 (d) x 1998 (*2076) (h) 1528(w) x 500 (d) x 1998 (*2076) (h) 1742 (w) x 500 (d) x 1998 (*2076) (h) 1314 (w) x 888 (d) x 1998 (*2081) (h) 1528 (w) x 888 (d) x 1998 (*2081) (h)   6(w) x10(h) 7(w) x10(h) 8(w) x10(h) 6(w) x10(h) 7(w) x10(h)   H-Temp, (Polysulfone) 3 years Warranty X-Temp, (Polyphenylsulfone) 8 years Warranty 1000 mm 1000 mm









## Quotation Ref: QUO-10219-R7C8X6

GM500 Components				
Product Code	Description	Unit Price	Qty	Total Price
GM500SU	GM500 H-Temp Polysulfone Cage base with releasing buttons & silicone gasket	£23.83	58	£1,382.14
GM500LID	GM500 Stainless Steel half wire bar Lid	£15.07	20	£301.40
GM500TOPSEALSU	GM500 H-Temp Polysulfone Cage Top with Bottle Depression, Microbiological Filter and Filter Retainer	£53.39	15	£800.85
ACBT0262SU	260 ml H-Temp Polysulfone Water Bottle with Silicone Ring	£4.20	220	£924.00
ACCP2521	S/S cap - 25 mm extension - hole Ø 1.8 mm	£5.77	130	£750.10
	Т	otal Price (e	excl. VA	AT) £4,158.49

	1284L Lids			
Product Code	Description	Unit Price	Qty	<b>Total Price</b>
1284L441SU	H-Temp Polysulfone Sealsafe Top with Bottle Depression, Filter and Clamps	£82.44	30	£2,473.20
		Total Price (e	xcl. V	AT) £2,473.20

Connections & Logistics				
Product Code	Description	Unit Price	Qty	Total Price
ACSCVF75M11UG	1.1m long 75mm diameter black autoclavable connection to connect the AHU to the IVC rack on the bottom left corner	£91.03	5	£455.15
9LBCBMP18SG3R	BC-Matic PLASTIC crate for 18 bottles 260ml square (TP ACBT0262) inclusive of grey plastic sliding lid	£221.87	4	£887.48
		Total Price (e	excl. VA	AT) £1,342.63

www.tecniplastuk.com



#### Page | 4

DGM racks, connections, AHU and Change Station				
Product Code	Description	Unit Price	Qty	Total Price
DGM70CPSU	Stock Clearance Single Sided DGM GM500 IVC Rack for 70 cages (7W X 10H) 1528 x 500 x 1998 (2076 with plenum) mm. Digital ready, Pre- assembled with 70 H-Temp Polysulfone GM500 cage bases, bottle filter tops and stainless steel half wire bar lids			
ACBT0262SU	<b>Stock Clearance</b> 260 ml H-Temp Polysulfone Water Bottle with Silicone Ring x 70 qty	Discounted Price £6,500.00	2	£13,000.00
ACCP3421GM	<b>Stock Clearance</b> Stainless Steel Bottle Cap with silicone seal x 70 qty			
ACPCP24GM	<b>Stock Clearance</b> Green Plastic Universal Card Holder x 70 qty			
ACSCVF33UG	Flexible "in line" rack interconnection kit	£127.62	4	£510.48
EMBOXUNEFUK	Universal Easy Flow for all IVC racks - 220/240V 50/60Hz (equipped with UK Plug)	£6,439.84	1	£6,439.84
9CS602TNNN	Aria CS60 larger open Change Station with reversible Trespa work surface 2 sockets, dimmable LED light	£10,479.57	1	£10,479.57

Estimated delivery time: 10 to 12 working weeks from receipt of purchase order.

### Free of charge delivery with this order

www.tecniplastuk.com



Page | 5

### **TECNIPLAST SPECIAL TERMS AND CONDITIONS**

- 1. Unless stated otherwise, payment terms are as follows: 100% payment for goods and services is due 30 days from the date of invoice or when available to deliver.
- 2. All prices are subject to VAT charged at the prevailing rate at time of order placement (valid exemption certificate accepted).
- 3. All prices are fixed and valid for 30 days.
- 4. The Stock Clearance items are sold as seen and carry no warranty.
- 5. The equipment is warranted against faulty parts and defects in workmanship for a period of 12 months from date of delivery or availability to deliver. The guarantee will not cover damage caused by accident or misuse through fault or negligence.
  - Warranty periods are:-
  - Racks 12 months.
  - AHU 24 months.
  - Cages (see point 6).
  - Laminar Flow 24 months.
- 6. In order for the manufacturer's warranty to remain valid for the full period, a routine maintenance must be performed by competent, qualified and fully trained personnel at the recommended intervals. Service Contracts can be quoted separately on request.

#### 7. Plastic Warranty:

• **TECNIPLAST H-TEMP Polysulfone (PSU)** caging is unconditionally guaranteed for **3** years from the stamped date.

During the warranty period Tecniplast liability shall be limited to the replacement of all severely warped, discoloured or cracked products not caused by improper use.

- Hereby follows a non-exhaustive list of improper usage examples:
- Washing with highly alkaline Ph detergents (particularly prolonged exposure to strong alkali in a hot environment).
- Autoclaving to temperatures higher than 134°C/273°F.
- Autoclaving in non-suitable plastic bags (sealed, no high steam permeability, wrong dimensions).

#### For full Terms and Conditions, please see our website <u>www.tecniplastuk.com</u>